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## ENGLISH

### Intended use:

APIS is designed to connect a syringe of ARCIOLANE (fractionated silicone oil for retina surgery) to silicone oil injection equipment.

### Presentation:

APIS is composed of: • 1 tube of 180 cm length • 1 white adaptor for the syringe • 1 male luer lock connection to compressed air source. • For equipment which have no male Luer compatible connection, an adapter is available according to the type of equipment:

Adapter reference	Equipment
No adapter	Premiere (Storz), Sto (Alcon), Penthesis (Ioltech), Antharès (Opticon)
A571	Associate, Harmony (Dorc) OpenPhaco (Corneal)
A572	Accurus (Alcon)
A573	Millenium (B&L)
A574	Orbit (Oërtli)
A575	Megatron S3 (Geuder)

The components are packaged in a NON sterilized labeled pouch. Instructions for use leaflet is attached. APIS can be re-sterilized with steam autoclave (see Sterilization). NOTE: adapters cannot be re-sterilized and can remain on the equipment between two operations.

### Instruction for use:

Only use the sterile APIS line. • Connect the syringe to the APIS connexion kit. • NB: Remove the white stopper block from the top of the syringe by rotating it a quarter turn • Insert the plug of the adaptor into the syringe barrel, • Adapt the ring on the upper part of the syringe flange. Check all connections are properly sealed. • Plug the other end of the APIS to a silicone oil injection machine by the means of the specific connector ( see above table ). • Injection of ARCIOLANE: • Remove the tip cap of the syringe and connect the desired surgical cannula to the syringe. Adjust the injection machine pressure.

### Precautions:

Maximum pressure allowed is 4 Bars. Do not handle the connectors when the line is under pressure.

### Sterilization :

Method of sterilization : autoclave 121° C, 20mm. The APIS line can be re-sterilized according to the following procedure : • At the end of each operation place the APIS line and its connection ring in demineralized water to moisten the inner tube. • Remove water excess by injecting air with a syringe through the cap hole. • Place the tube and the connection ring into a steam sterilization compatible pouch-bag. • Seal the pouch bag. • Put the pouched product in steam autoclave. • Select a cycle with a plateau at 121° C for 20 minutes. • When cycle is over , make sure the product is dry before removing from autoclave and let the product return to ambient temperature. • Keep APIS away from aggressive chemical compounds such as sodium hydroxide, hydrochloric acid, bleach, formaldehyde...

### Storage conditions:

At room temperature and in a dry place.

Manufactured by ARCADOPHTA